

**ANSI Tipo K 11.8kg Pull Power (herradura)  
Termopar magnético, cable aislado PFA con trenza  
de acero inoxidable que termina en miniatura o  
conchete estándar**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Características técnicas****Specifications**

<b>General Description</b>	11.8kg Pull Power (Horseshoe) Magnet Thermocouple
<b>Sensor Type</b>	Magnet Thermocouple
<b>Thermocouple Type</b>	K
<b>Magnet</b>	45 x 30 x 30mm (L x H x W) - Max pull: 11.8 KG
<b>Cable Length</b>	1.5m or 3m
<b>Cable Type</b>	PFA Insulated Cable with Stainless Steel Over-Braid
<b>Core / Strands</b>	7/0.2mm
<b>Cable Termination</b>	Miniature or Standard Plug
<b>Max. Temperature</b>	+250°C
<b>Min. Temperature</b>	-50°C
<b>Standards Met</b>	ANSI MC96.1