



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Termopar magnético de potencia de tracción (herradura) tipo K JIS, cable aislado de PFA con trenza de acero inoxidable que termina en



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Caracterstias tcnicas**Specifications**

Product Code	XE-0674-001
General Description	4.5kg Pull Power (Horseshoe) Magnet Thermocouple
Sensor Type	K JIS
Magnet	30 x 20 x 20mm (L x H x W). Max Pull: 4.5kg
Cable Length	1M
Cable Type	Glassfibre insulated lead with stainless steel overbraid
Core / Strands	7/0.2mm
Termination Type	Miniature Plug
Max. Temperature	+250°C
Min. Temperature	-50°C
Standards Met	JIS (Japanese Industrial Standards)