


**Termopar magnético de potencia de tracción  
(herradura) tipo K JIS de 4,5 kg, cable aislado con  
PFA con trenza de acero inoxidable que termina  
en enchufe estándar - 2 m**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

## Características técnicas

**Specifications**

<b>Product Code</b>	XE-0677-001
<b>General Description</b>	4.5kg Pull Power (Horseshoe) Magnet Thermocouple
<b>Sensor Type</b>	K JIS
<b>Magnet</b>	30 x 20 x 20mm (L x H x W). Max Pull: 4.5kg
<b>Cable Length</b>	2m
<b>Cable Type</b>	Glassfibre insulated lead with stainless steel overbraid
<b>Core / Strands</b>	7/0.2mm
<b>Termination Type</b>	Standard Plug
<b>Max. Temperature</b>	250°C
<b>Min. Temperature</b>	-50°C
<b>Standards Met</b>	JIS (Japanese Industrial Standards)