



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

**JIS Tipo K 4.5kg Potencia de traccin (herradura)
Termopar magnitico, cable aislado PFA con trenza
de acero inoxidable que termina en miniatura o
conchife estándar**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Características técnicas**Specifications**

General Description	4.5kg Pull Power (Horseshoe) Magnet Thermocouple
Sensor Type	Magnet Thermocouple
Thermocouple Type	K
Magnet	30 x 20 x 20mm (L x H x W). Max Pull: 4.5kg
Cable Length	1.0m or 2.0m
Cable Type	Glassfibre insulated lead with stainless steel overbraid
Core / Strands	7/0.2mm
Termination Type	Miniature or Standard Plug
Max. Temperature	+250°C
Min. Temperature	-50°C
Standards Met	JIS (Japanese Industrial Standards)