

**Termopar magnético de potencia de traccin  
(herradura) tipo K JIS, cable aislado de PFA con  
trenza de acero inoxidable que termina en  
enchufe en miniatura - 3 m**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Características técnicas****Specifications**

<b>Product Code</b>	XE-0679-001
<b>General Description</b>	Rapid, reliable temperature sensing on ferrous surfaces
<b>Sensor Type</b>	K JIS
<b>Thermocouple Junction</b>	Grounded at Tip
<b>Magnet</b>	25 x 25 x 40mm (HxWxL) with powerful 9kg pull
<b>Cable Length</b>	3m
<b>Cable Type</b>	PFA Teflon® insulated with stainless steel over-braid, fitted plug termination
<b>Termination Type</b>	Miniature Plug
<b>Max. Temperature</b>	250°C
<b>Min. Temperature</b>	-50°C
<b>Standards Met</b>	JIS (Japanese Industrial Standards)