


**Termopar magnético de potencia de tracción  
(herradura) tipo K JIS Tipo K, cable aislado con  
PFA con trenza de acero inoxidable que termina  
en enchufe en miniatura - 1,5 m**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Caracterstias tcnicas****Specifications**

<b>Product Code</b>	XE-0682-001
<b>General Description</b>	11.8kg Pull Power (Horseshoe) Magnet Thermocouple
<b>Sensor Type</b>	Magnet Thermocouple
<b>Thermocouple Type</b>	K JIS
<b>Magnet</b>	45 x 30 x 30mm (L x H x W) - Max pull: 11.8 KG
<b>Cable Length</b>	1.5m
<b>Cable Type</b>	PFA Insulated Cable with Stainless Steel Over-Braid
<b>Core / Strands</b>	7/0.2mm
<b>Termination Type</b>	Miniature Plug
<b>Max. Temperature</b>	250°C
<b>Min. Temperature</b>	-50°C
<b>Standards Met</b>	JIS (Japanese Industrial Standards)