

Termopar magnético de potencia de tracción (herradura) tipo K JIS de 11,8 kg, cable aislado de PFA con trenza de acero inoxidable que termina en enchufe en miniatura 2 m



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Caracterstias tcnicas**Specifications**

Product Code	XE-0683-001
General Description	11.8kg Pull Power (Horseshoe) Magnet Thermocouple
Sensor Type	Magnet Thermocouple
Thermocouple Type	K JIS
Magnet	45 x 30 x 30mm (L x H x W) - Max pull: 11.8 KG
Cable Length	3m
Cable Type	PFA Insulated Cable with Stainless Steel Over-Braid
Core / Strands	7/0.2mm
Termination Type	Miniature Plug
Max. Temperature	250°C
Min. Temperature	-50°C
Standards Met	JIS (Japanese Industrial Standards)