



Conector de fascia de termopar rectangular estndar IS-J-FF Tipo J IEC



Standard Fascia Socket (with nylon mounting clip)
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Conector de fascia de termopar rectangular estndar IS-J-FF Tipo J IEC

Clip de nylon para una fijacin rpida y sencilla al panel. Todos los zcalos de montaje de la fascia estn polarizados por medio de pasadores de tamao desigual

Caractersticas tcnicas**Specifications**

| | |
|----------------------------|---|
| Product Code | XE-1467-001 |
| General Description | Accepts standard round pin plug connector. Nylon clip for quick, simple fixing to panel |
| Thermocouple Type | J IEC |
| Connector Type | Fascia Mounting Socket |
| Connector Size | Standard |
| Colour | Black |
| Max. Temperature | 120°C |
| +Positive Contact | Iron |
| -Negative Contact | Copper Nickel |
| Standards Met | BS EN 50212:1997. RoHS |
| Series | Standard Connector |