



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Termopar de hoja JIS, cable aislado en fibra de vidrio con trenza de acero inoxidable que termina en colas desnudas, enchufe en miniatura o



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Caractersticas tcnicas**Specifications**

General Description	Thermocouple mounted on a stainless steel shim
Thermocouple Type	K
Thermocouple Junction	Junction grounded at tip for fast response
Cable Length	1,2,3 or 5m
Cable Type	Glassfibre insulated lead with stainless steel overbraid
Core / Strands	7/0.2mm
Termination Type	Bare tails, Miniature or Standard Plug (with cable clamps)
Max. Temperature	+350°C
Min. Temperature	-60°C
Standards Met	JIS (Japanese Industrial Standards)
Dimensions	25x13mm shim (0.2mm thick)