



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Sonda de termopar de uso general ANSI, cable aislado en fibra de vidrio con trenza de acero inoxidable que termina en colas desnudas, enchufe en miniatura a estándar Tipo K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Caractersticas tcnicas**Specifications**

General Description	Rugged thermocouple with stainless steel sheath
Thermocouple Type	K
Thermocouple Junction	Grounded hot junction for fast response
Probe Diameter	4.5mm
Probe Length	150mm
Cable Length	1 or 2m
Cable Type	Glassfibre insulated lead with stainless steel overbraid
Core / Strands	7/0.2mm
Termination Type	Bare tails, Miniature or Standard Plug (with cable clamps)
Max. Temperature	+350°C
Min. Temperature	-60°C
Standards Met	ANSI MC96.1