


**Sonda de termopar de uso general, cable aislado en fibra de vidrio con trenza de acero inoxidable - Tipo K JIS con terminacin de enchufe estndar - 2 m de longitud**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Caractersticas tcnicas****Specifications**

<b>Product Code</b>	XE-0771-001
<b>General Description</b>	Rugged thermocouple with stainless steel sheath
<b>Thermocouple Type</b>	K JIS
<b>Thermocouple Junction</b>	Grounded hot junction for fast response
<b>Probe Diameter</b>	4.5mm
<b>Probe Length</b>	150mm
<b>Cable Length</b>	2M
<b>Cable Type</b>	Glassfibre insulated lead with stainless steel overbraid
<b>Core / Strands</b>	7/0.2mm
<b>Termination Type</b>	Standard Plug (with cable clamps)
<b>Max. Temperature</b>	+350°C
<b>Min. Temperature</b>	-60°C
<b>Standards Met</b>	JIS (Japanese Industrial Standards)