

**Termopar de boquilla ANSI, cable aislado en fibra de vidrio con trenza de acero inoxidable que termina en colas desnudas, enchufe en miniatura**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Caractersticas tcnicas****Specifications**

<b>General Description</b>	Thermocouple housed in a stainless steel threaded bolt
<b>Thermocouple Type</b>	K
<b>Thermocouple Junction</b>	Grounded hot junction for fast response, M6x1mm thread
<b>Cable Length</b>	1, 2 or 5m
<b>Cable Type</b>	Glassfibre insulated lead with stainless steel overbraid
<b>Core / Strands</b>	7/0.2mm
<b>Termination Type</b>	Bare tails, Miniature or Standard Plug (with cable clamps)
<b>Max. Temperature</b>	+350°C
<b>Min. Temperature</b>	-60°C
<b>Standards Met</b>	ANSI MC96.1