

**Termopar de boquilla, cable aislado en fibra de vidrio con trenza de acero inoxidable - Tipo K JIS con terminacin de cola desnuda - 1 m de longitud**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Caracterstias tcnicas****Specifications**

<b>Product Code</b>	XE-0893-001
<b>General Description</b>	Thermocouple housed in a stainless steel threaded bolt
<b>Thermocouple Type</b>	K JIS
<b>Thermocouple Junction</b>	Grounded hot junction for fast response, M6x1mm thread
<b>Cable Length</b>	1M
<b>Cable Type</b>	Glassfibre insulated lead with stainless steel overbraid
<b>Core / Strands</b>	7/0.2mm
<b>Termination Type</b>	Bare Tails
<b>Max. Temperature</b>	+350°C
<b>Min. Temperature</b>	-60°C
<b>Standards Met</b>	JIS (Japanese Industrial Standards)