

**Termopar de arandela ANSI, cable aislado en fibra de vidrio con trenza de acero inoxidable que termina en colas desnudas o enchufe estndar -**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Caractersticas tcnicas****Specifications**

<b>General Description</b>	Thermocouple housed in a washer eyelet
<b>Thermocouple Type</b>	K
<b>Thermocouple Junction</b>	Grounded tip for quick response
<b>Cable Length</b>	2m
<b>Cable Type</b>	Glassfibre insulated, stainless steel overbraided
<b>Core / Strands</b>	7/0.2mm
<b>Termination Type</b>	Either tails or standard plug termination with cable clamps
<b>Max. Temperature</b>	+350°C
<b>Standards Met</b>	ANSI
<b>Washer</b>	3.5mm, 4.0mm, 5.0mm, 6.0mm, 8.0mm & 10.0mm clearance hole available