



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

**Termopar de arandela, cable aislado de fibra de vidrio con trenza de acero inoxidable - Tipo K JIS con terminacin de enchufe estndar y tamao de orificio de arandela de 4,0 mm**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Caractersticas tcnicas****Specifications**

<b>Product Code</b>	XE-0972-001
<b>General Description</b>	Thermocouple housed in a washer eyelet
<b>Thermocouple Type</b>	K JIS
<b>Thermocouple Junction</b>	Grounded tip for quick response
<b>Cable Length</b>	2m
<b>Cable Type</b>	Glassfibre insulated, stainless steel overbraided
<b>Core / Strands</b>	7/02.mm
<b>Termination Type</b>	Standard plug termination with cable clamps
<b>Max. Temperature</b>	+350°C
<b>Standards Met</b>	JIS (Japanese Industrial Standards)
<b>Washer</b>	4.0mm clearance hole