



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Termopar de arandela JIS, cable aislado en fibra de vidrio con trenza de acero inoxidable que termina en colas desnudas o enchufe estndar -



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Caractersticas tcnicas**Specifications**

General Description	Thermocouple housed in a washer eyelet
Thermocouple Type	K
Thermocouple Junction	Grounded tip for quick response
Cable Length	2m
Cable Type	Glassfibre insulated, stainless steel overbraided
Core / Strands	7/02.mm
Termination Type	Either tails or standard plug termination with cable clamps
Max. Temperature	+350°C
Standards Met	JIS
Washer	3.5mm, 4.0mm, 5.0mm, 6.0mm, 8.0mm & 10.0mm clearance hole available