

Aceitunas PTFE



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Aceitunas PTFE. Disponible para adaptarse a dimetros de sonda que van desde 1,5 mm a 8,0 mm

Specifications

General Description	Olives
Material	PTFE
Probe Diameter	Suitable for probe diameters between 1.5mm & 8.0mm