

Sonda de punta plana - Tipo RTD / PRT - 2M de longitud



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Sonda Pt100 de uso general

Aplanado en la punta para mejorar la respuesta trmica

Configuracin de 2 hilos

Clase B, segn IEC 751

Funda de acero inoxidable 316, 4,0 mm de dimetro x 150 mm de largo

Cable trenzado de 2 ncleos con aislamiento de tefln de 7/0,2 mm de 1 metro de largo

Incluye accesorio de compresin macho de latn BSPP de 1/8" con tuerca de bloqueo

Rango de temperatura de -50 C a + 250 C

Caractersticas tcnicas

Specifications

Product Code	XE-3698-001
General Description	316 stainless steel probe, 4.0mm diameter (6mm max. across flat area of tip) x 150mm long overall
Sensor Type	Pt100 (100 Ohms @ 0°C) to IEC 751, Class B
Probe Temperature Range	-50°C to +250°C
Cable Length	2M
Cable Type	7/0.2mm Teflon® insulated 2 core twisted leads
Termination Type	2-wire tails
Accessories	1/8" BSPP brass male compression fitting with locknut included