

Conector de termopar en miniatura Zcalo en Inea BM-K-FS Tipo K BS



Miniature in-line Socket
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Conector de termopar en miniatura Zcalo en Inea BM-K-FS Tipo K BS

Ideal para uso general, todos los contactos estn polarizados para garantizar una conexin correcta. Estos conectores aceptarn cables de termopar de hasta 4 mm de dimetro.

Caractersticas tcnicas**Specifications**

Product Code	XE-1386-001
General Description	Ideal for general purpose use. Contacts are polarised to ensure correct connection.
Thermocouple Type	K BS
Connector Type	Socket
Connector Size	Miniature
Colour	Yellow
Max. Temperature	220°C
+Positive Contact	Nickel Chromium
-Negative Contact	Nickel Alloy
Standards Met	BS EN 50212:1997. RoHS