


## Bloque de terminal de barrera Conector de termopar IS-CU-BTS Tipo CU



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Bloque de terminal de barrera Conector de termopar IS-CU-BTS Tipo CU**

Se suministra en pares individuales entrelazados. Conexiones en material termopar

**Especificaciones****Specifications**

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | XE-1944-001  |
| <b>General Description</b> | Supplied in single, interlocking pairs. Connections in thermocouple material |
| <b>Thermocouple Type</b>   | IEC, ANSI, JIS Type CU   |
| <b>Connector Type</b>      | Barrier Terminal Blocks  |
| <b>Colour</b>              | White  |
| <b>Dimensions</b>          | See Schematic Drawing (Product Literature Tab)                               |
| <b>Max. Temperature</b>    | 220°C  |
| <b>+Positive Contact</b>   | Copper   |
| <b>-Negative Contact</b>   | Copper   |
| <b>Standards Met</b>       | BS EN 50212:1997. RoHS   |