


Bloque de terminal de barrera Conector de termopar IS-CU-BTS Tipo CU



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Bloque de terminal de barrera Conector de termopar IS-CU-BTS Tipo CU

Se suministra en pares individuales entrelazados. Conexiones en material termopar

Especificaciones**Specifications**

Product Code	XE-1944-001
General Description	Supplied in single, interlocking pairs. Connections in thermocouple material
Thermocouple Type	IEC, ANSI, JIS Type CU
Connector Type	Barrier Terminal Blocks
Colour	White
Dimensions	See Schematic Drawing (Product Literature Tab)
Max. Temperature	220°C
+Positive Contact	Copper
-Negative Contact	Copper
Standards Met	BS EN 50212:1997. RoHS