



Terminales de pala de termopar tipo K



Thermocouple Spade Terminal

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Terminales de pala engarzadosPrecio por paquete de 20 paresAdecuado para los tipos de calibracin IEC, ANSI y JIS.Fabricado con material de termopar compatibleTipo de termopar estampado en terminalSuits Labfacility & otros terminales de barreraPara tamaos de cable de 0,5 a 1,0 mm<p></p>

Specifications

Product Code	XE-1982-020
General Description	pack of 20 pairs
+Positive Contact	Nickel Chromium
-Negative Contact	Nickel Alloy