



## Terminales de pala de termopar tipo K



Thermocouple Spade Terminal

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

<ul><li>Terminales de pala engarzados</li><li>Precio por paquete de 20 pares</li><li>Adecuado para los tipos de calibracin IEC, ANSI y JIS.</li><li>Fabricado con material de termopar compatible</li><li>Tipo de termopar estampado en terminal</li><li>Suits Labfacility & otros terminales de barrera</li><li>Para tamaos de cable de 0,5 a 1,0 mm</li></ul><p></p>

## Specifications

|                            |                  |
|----------------------------|------------------|
| <b>Product Code</b>        | XE-1982-020      |
| <b>General Description</b> | pack of 20 pairs |
| <b>+Positive Contact</b>   | Nickel Chromium  |
| <b>-Negative Contact</b>   | Nickel Alloy     |