




## Terminales de pala termopar tipo CU



Thermocouple Spade Terminal

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Terminales de pala engarzados

Precio por pack de 20 pares

Adecuado para los tipos de calibración IEC, ANSI y JIS.

Hecho de material de termopar compatible

Tipo de termopar estampado en el terminal

Se adapta a Labfacility y otros terminales de barrera

Para tamanos de alambre de 0,5 a 1,0 mm+contacto positivo

## Specifications

<b>Product Code</b>	XE-1986-020
<b>General Description</b>	pack of 20 pairs
<b>+Positive Contact</b>	Copper
<b>-Negative Contact</b>	Copper