


Conector de termopar en miniatura Am-N-M Tipo N ANSI



Miniature in-line Plug
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Conector de termopar en miniatura Am-N-M Tipo N ANSI

Con pasadores planos para adaptarse al tipo de termopar e ideal para uso general, todos los contactos están polarizados para garantizar una correcta conexión. Estos conectores aceptarán un cable termopar de hasta 4 mm de diámetro.

Especificaciones**Specifications**

Product Code	XE-1270-001
General Description	Ideal for general purpose use. Contacts are polarised to ensure correct connection.
Thermocouple Type	N ANSI
Connector Type	Plug
Connector Size	Miniature
Colour	Orange
Max. Temperature	220°C
Pins	Flat Pins
+Positive Contact	Nickel / Chromium / Silicon
-Negative Contact	Nickel Silicon
Standards Met	BS EN 50212:1997. RoHS