

## Conector de termopar en miniatura Zcalo en Inea AM-K-FS Tipo K ANSI



**Miniature in-line Socket**  
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Conector de termopar en miniatura Zcalo en lnea AM-K-FS Tipo K ANSI**

Ideal para uso general, todos los contactos estn polarizados para garantizar una conexin correcta. Estos conectores aceptarn cables de termopar de hasta 4 mm de dimetro.

**Caractersticas tcnicas****Specifications**

<b>Product Code</b>	XE-1132-001
<b>General Description</b>	Ideal for general purpose use. Contacts are polarised to ensure correct connection.
<b>Thermocouple Type</b>	K ANSI
<b>Connector Type</b>	Socket
<b>Connector Size</b>	Miniature
<b>Colour</b>	Yellow
<b>Max. Temperature</b>	220°C
<b>+Positive Contact</b>	Nickel Chromium
<b>-Negative Contact</b>	Nickel Alloy
<b>Standards Met</b>	BS EN 50212:1997. RoHS