



Conector de termopar en miniatura Zcalo en Inea AM-T-FS Tipo T ANSI



Miniature in-line Socket
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Conector de termopar en miniatura Zcalo en Inea AM-T-FS Tipo T ANSI

Ideal para uso general, todos los contactos estn polarizados para garantizar una conexin correcta. Estos conectores aceptarn cables de termopar de hasta 4 mm de dimetro.

Caractersticas tcnicas**Specifications**

Product Code	XE-1222-001
General Description	Ideal for general purpose use. Contacts are polarised to ensure correct connection.
Thermocouple Type	T ANSI
Connector Type	Socket
Connector Size	Miniature
Colour	Blue
Max. Temperature	220°C
+Positive Contact	Copper
-Negative Contact	Copper Nickel
Standards Met	BS EN 50212:1997. RoHS
Series	Miniature Connector