


## Conector de termopar rápido en miniatura AM-K-MQ Tipo K ANSI



**Miniature in-line Plug - Quick Wire**  
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Conector de termopar en miniatura Quick Wire Plug AM-K-MQ Tipo K ANSI**

Conector en miniatura con pasadores de enchufe planos. Fácil conexión de jab-in para una terminación rápida. Sólo empuje el cable y apriete el tornillo.

**Especificaciones****Specifications**

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | XE-1141-001                                   |
| <b>General Description</b> | Easy jab-in connection for rapid termination. |
| <b>Thermocouple Type</b>   | K ANSI  |
| <b>Connector Type</b>      | Quick Wire Plug                               |
| <b>Connector Size</b>      | Miniature                                     |
| <b>Colour</b>              | Yellow  |
| <b>Max. Temperature</b>    | 220°C   |
| <b>Pins</b>                | Flat Pins                                     |
| <b>+Positive Contact</b>   | Nickel Chromium                               |
| <b>-Negative Contact</b>   | Nickel Alloy                                  |
| <b>Standards Met</b>       | BS EN 50212:1997. RoHS                        |