



Conector de termopar rápido en miniatura AM-N-MQ Tipo N ANSI



Miniature in-line Plug - Quick Wire
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Conector de termopar en miniatura Quick Wire Plug AM-N-MQ Tipo N ANSI

Conector en miniatura con pasadores de enchufe planos. Fácil conexión de jab-in para una terminación rápida. Sólo empuje el cable y apriete el tornillo.

Especificaciones**Specifications**

Product Code	XE-1281-001
General Description	Easy jab-in connection for rapid termination.
Thermocouple Type	N ANSI
Connector Type	Quick Wire Plug
Connector Size	Miniature
Colour	Orange
Max. Temperature	220°C
Pins	Flat Pins
+Positive Contact	Nickel / Chromium / Silicon
-Negative Contact	Nickel Silicon
Standards Met	BS EN 50212:1997. RoHS