

## WT9055 Sonmetro digital



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

El sonometro digital WT9055 es novedoso, pequeno y porttil. El dispositivo es aplicable para la medicin de entornos de fbrica, control de calidad, prevencin de la salud y diversos ruidos ambientales, incluidas oficinas, rutas de transporte, equipos estreo y muchos ms.

Medicin del nivel sonoro (dBA)

MIN / MAX / Valor actual de bloqueo

Retencin de datos

Funcin de retroiluminacin LCD

Apagado manual / automtico

Alarma de retroiluminacin

### Caractersticas tcnicas

## Specifications

<b>Product Code</b>	XE-4157-001
<b>Model</b>	WT9055
<b>Type</b>	Sound Level Meter
<b>Accuracy</b>	+/- 1.5dB
<b>Resolution</b>	0.1dB(A)
<b>Battery</b>	3 x 1.5V AAA
<b>Dimensions</b>	53.4mm x 31.4mm x 131.7mm
<b>Working Temperature</b>	0-40°C
<b>Operating Humidity</b>	10-80% RH
<b>Weight</b>	96.38g
<b>Storage Temperature</b>	10-60°C
<b>Measuring range</b>	35-130dB (A)
<b>Important Information</b>	Frequency response: 31.5Hz - 8KHz