



Conector de termopar de zcalo de fascia rectangular en miniatura AM-K-FF Tipo K ANSI



Miniature Fascia Socket (with nylon mounting clip)
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Conector de termopar rectangular en miniatura Zcalo de fascia AM-K-FF Tipo K ANSI

Clip de nylon para una fijacin rpida y sencilla al panel. Todos los zcalos de montaje de la fascia estn polarizados por medio de clavijas de tamao desigual

Caractersticas tcnicas**Specifications**

Product Code	XE-1137-001
General Description	All fascia mounting sockets are polarised by means of unequally sized pins
Thermocouple Type	K ANSI
Connector Type	Fascia mounting socket
Connector Size	Miniature
Colour	Yellow
Dimensions	See Schematic Drawing (Product Literature Tab)
Max. Temperature	120°C
+Positive Contact	Nickel Chromium
-Negative Contact	Nickel Alloy
Standards Met	BS EN 50212:1997. RoHS