

Miniatura Rectangular Fascia Socket Conector termopar AM-R/S-FF Tipo R/S ANSI



Miniatura Fascia Socket (with nylon mounting clip)
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Conector de termopar rectangular en miniatura Fascia Socket AM-R/S-FF Tipo R/S ANSI

Clip de nylon para una fijación rápida y sencilla al panel. Todas las tomas de montaje de fascia están polarizadas por medio de pasadores de igual tamaño

Especificaciones**Specifications**

Product Code	XE-1317-001
General Description	All fascia mounting sockets are polarised by means of unequally sized pins
Thermocouple Type	R/S ANSI
Connector Type	Fascia mounting socket
Connector Size	Miniature
Colour	Green
Max. Temperature	120°C
+Positive Contact	Copper
-Negative Contact	Copper Alloy
Standards Met	BS EN 50212:1997. RoHS