

## Miniatura Rectangular Fascia Socket Conector termopar AM-E-FF Tipo E ANSI



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Conector de termopar rectangular en miniatura Fascia Socket AM-E-FF Tipo E ANSI**

Clip de nylon para una fijación rápida y sencilla al panel. Todas las tomas de montaje de fascia están polarizadas por medio de pasadores de igual tamaño

**Especificaciones****Specifications**

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | XE-9037-001  |
| <b>General Description</b> | All fascia mounting sockets are polarised by means of unequally sized pins |
| <b>Thermocouple Type</b>   | E ANSI   |
| <b>Connector Type</b>      | Fascia mounting socket   |
| <b>Connector Size</b>      | Miniature  |
| <b>Colour</b>              | Violet   |
| <b>Dimensions</b>          | See Schematic Drawing (Product Literature Tab)                             |
| <b>Max. Temperature</b>    | 120°C  |
| <b>+Positive Contact</b>   | Nickel Chromium  |
| <b>-Negative Contact</b>   | Copper Nickel  |
| <b>Standards Met</b>       | BS EN 50212:1997. RoHS   |