



Conector de termopar de montaje de PCB Socket AM-T-PCB Tipo T ANSI



Identification sticker to denote calibration type: K, J, T, N, E, R/S or CU



Miniature PCB Socket
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Conector de termopar de montaje de PCB Socket AM-T-PCB Tipo T ANSI

Especialmente diseñado para PCBs, pines sólidos para montaje y aprovisionamiento de placa directa para sensor CJC

Especificaciones**Specifications**

Product Code	XE-1229-001
General Description	Versatile, high quality circuit board mounting miniature socket
Thermocouple Type	T ANSI
Connector Type	PCB socket
Connector Size	Miniature
Colour	Blue
Max. Temperature	220°C
+Positive Contact	Copper
-Negative Contact	Copper Nickel
Standards Met	BS EN 50212:1997. RoHS
Series	Miniature Connector