



Conector de termopar de montaje de PCB Socket AM-T-PCB Tipo T ANSI



Miniature PCB Socket
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Conector de termopar de montaje de PCB Socket AM-T-PCB Tipo T ANSI

Especialmente diseñado para PCBs, pines sólidos para montaje y aprovisionamiento de placa directa para sensor CJC

Especificaciones**Specifications**

| | |
|----------------------------|---|
| Product Code | XE-1229-001 |
| General Description | Versatile, high quality circuit board mounting miniature socket |
| Thermocouple Type | T ANSI |
| Connector Type | PCB socket |
| Connector Size | Miniature |
| Colour | Blue |
| Max. Temperature | 220°C |
| +Positive Contact | Copper |
| -Negative Contact | Copper Nickel |
| Standards Met | BS EN 50212:1997. RoHS |