

## Enchufe de conector de termopar estndar AS-K-M Tipo K ANSI



**Standard in-line Plug**  
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Enchufe de conector de termopar estndar AS-K-M Tipo K ANSI**

Los contactos estn hechos de verdaderas aleaciones de termopar. Adecuado para aplicaciones industriales

**Caractersticas tcnicas****Specifications**

<b>Product Code</b>	XE-1430-001
<b>General Description</b>	Contacts are polarised to prevent incorrect connection
<b>Thermocouple Type</b>	K ANSI
<b>Connector Type</b>	Plug
<b>Connector Size</b>	Standard
<b>Colour</b>	Yellow
<b>Dimensions</b>	See Schematic Drawing (Product Literature Tab)
<b>Max. Temperature</b>	220°C
<b>Pins</b>	Solid Round Pins
<b>+Positive Contact</b>	Nickel Chromium
<b>-Negative Contact</b>	Nickel Alloy
<b>Standards Met</b>	BS EN 50212:1997. RoHS