

# Conector de termopar estãfÂndar Conector en IAfÂ-nea As-E-F Tipo E ANSI



**Standard in-line Socket**  
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Conector de termopar estándar Conector en línea As-E-F Tipo E ANSI**

Acepta pasadores de enchufe redondos estándar, adecuados para aplicaciones industriales

**Especificaciones****Specifications**

<b>Product Code</b>	XE-9112-001
<b>General Description</b>	Accepts standard round plug pins, contacts are polarised to prevent incorrect connection.
<b>Thermocouple Type</b>	E ANSI
<b>Connector Type</b>	Socket
<b>Connector Size</b>	Standard
<b>Colour</b>	Violet
<b>Max. Temperature</b>	220°C
<b>+Positive Contact</b>	Nickel Chromium
<b>-Negative Contact</b>	Copper Nickel
<b>Standards Met</b>	BS EN 50212:1997. RoHS