


## Zcalo de conector de termopar de fascia rectangular estndar AS-K-FF Tipo K ANSI



**Standard Fascia Socket (with nylon mounting clip)**  
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Conector de termopar rectangular estndar Zcalo de fascia AS-K-FF Tipo K ANSI**

Pinza de nylon para una fijacin rpida y sencilla al panel. Todos los zcalos de montaje de la fascia estn polarizados por medio de clavijas de tamao desigual

**Caractersticas tcnicas****Specifications**

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | XE-1437-001   |
| <b>General Description</b> | Accepts standard round pin plug connector. Nylon clip for quick, simple fixing to panel |
| <b>Thermocouple Type</b>   | K ANSI  |
| <b>Connector Type</b>      | Fascia Mounting Socket  |
| <b>Connector Size</b>      | Standard  |
| <b>Colour</b>              | Yellow  |
| <b>Max. Temperature</b>    | 120°C   |
| <b>+Positive Contact</b>   | Nickel Chromium   |
| <b>-Negative Contact</b>   | Nickel Alloy  |
| <b>Standards Met</b>       | BS EN 50212:1997. RoHS  |