



Bloque de terminales de barrera Conector de termopar AS-T-BTS Tipo T ANSI



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Bloque de terminales de barrera Conector de termopar AS-T-BTS Tipo T ANSI

Se suministra en pares individuales entrelazados. Conexiones en material termopar

Especificaciones**Specifications**

Product Code	XE-1952-001
General Description	Supplied in single, interlocking pairs. Connections in thermocouple material
Thermocouple Type	T ANSI
Connector Type	Barrier Terminal Blocks
Colour	Blue
Dimensions	See Schematic Drawing (Product Literature Tab)
Max. Temperature	220°C
+Positive Contact	Copper
-Negative Contact	Copper Nickel
Standards Met	BS EN 50212:1997. RoHS