



# Conector de termopar de montaje pcb Conector JM-B-PCB Tipo B JIS



**Miniature PCB Socket**  
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Conector de termopar de montaje pcb Conector JM-B-PCB Tipo B JIS**

Especialmente diseñado para PCBs, pines sólidos para montaje y aprovisionamiento de placa directa para sensor CJC

**Especificaciones****Specifications**

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | XE-9220-001   |
| <b>General Description</b> | Versatile, high quality circuit board mounting miniature socket |
| <b>Thermocouple Type</b>   | B JIS   |
| <b>Connector Type</b>      | PCB Socket  |
| <b>Connector Size</b>      | Miniature   |
| <b>Colour</b>              | Grey  |
| <b>Max. Temperature</b>    | 220°C   |
| <b>+Positive Contact</b>   | Copper  |
| <b>-Negative Contact</b>   | Copper  |
| <b>Standards Met</b>       | BS EN 50212:1997. RoHS  |