



## Conector de termopar de montaje en panel estándar con soporte de acero inoxidable JS-E-SSPF Tipo E JIS



**Standard in-line Socket with stainless steel mounting bracket**  
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Conector de termopar de montaje en panel estándar con soporte de acero inoxidable JS-E-SSPF Tipo E JIS**

Destinado para su uso en procesos más industriales, pero se puede utilizar en todas las aplicaciones

**Especificaciones****Specifications**

<b>Product Code</b>	XE-9358-001
<b>General Description</b>	Intended for use in more industrial processes but can be used in all applications as desired
<b>Thermocouple Type</b>	E JIS
<b>Connector Type</b>	Socket with SS Bracket
<b>Connector Size</b>	Standard
<b>Colour</b>	Violet
<b>Max. Temperature</b>	220°C
<b>+Positive Contact</b>	Nickel Chromium
<b>-Negative Contact</b>	Copper Nickel
<b>Standards Met</b>	BS EN 50212:1997. RoHS