



## Bloque de terminales de barrera Conector de termopar JS-K-BTS Tipo K JIS



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Bloque de terminales de barrera Conector de termopar JS-K-BTS Tipo K JIS**

Se suministra en pares individuales entrelazados. Conexiones en material termopar

**Especificaciones****Specifications**

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | XE-1960-001  |
| <b>General Description</b> | Supplied in single, interlocking pairs. Connections in thermocouple material |
| <b>Thermocouple Type</b>   | K JIS  |
| <b>Connector Type</b>      | Barrier Terminal Blocks  |
| <b>Colour</b>              | Blue   |
| <b>Dimensions</b>          | See Schematic Drawing (Product Literature Tab)                               |
| <b>Max. Temperature</b>    | 220°C  |
| <b>+Positive Contact</b>   | Nickel Chromium  |
| <b>-Negative Contact</b>   | Nickel Alloy   |
| <b>Standards Met</b>       | BS EN 50212:1997. RoHS   |