


Ojal de alivio de tensin



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Arandela de alivio de tensin para conectores en miniatura

Rango de temperatura: -30 C a +140 C

Material: Elastmero termoplstico (TPE)

Specifications

| | |
|---------------------|-------------|
| Product Code | XE-1882-001 |
|---------------------|-------------|