



**LABFACILITY**

TEMPERATURE & PROCESS TECHNOLOGY

## Cable de termopar curlado retrctil tipo T (IEC)



XE-7033-001  
**TYPE-T**  
(IEC)



XE-7033-001  
**TYPE-T** (IEC)



XE-7033-001  
**TYPE-T** (IEC)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Disponible en tipo T

Codificado por colores segn IEC

2 metros de longitud cuando est completamente extendido

Terminado con un enchufe en miniatura instalado

Exterior de poliuretano

## Specifications

<b>Product Code</b>	XE-7033-001
<b>Cable Length</b>	2M fully extended
<b>Cable Type</b>	Polyurethane outer
<b>Connector Type</b>	Thermocouple Plug
<b>Connector Size</b>	Miniature
<b>Colour</b>	Brown
<b>Standards Met</b>	Colour code (cores & jacket) to IEC-584-3
<b>Thermocouple Type</b>	T
<b>Termination Type</b>	Miniature Thermocouple Plug