

Cable de termopar curlado retrctil tipo T (IEC)



XE-7033-001
TYPE-T
(IEC)



XE-7033-001
TYPE-T (IEC)



XE-7033-001
TYPE-T (IEC)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Disponible en tipo T

Codificado por colores segn IEC

2 metros de longitud cuando est completamente extendido

Terminado con un enchufe en miniatura instalado

Exterior de poliuretano

Specifications

Product Code	XE-7033-001
Cable Length	2M fully extended
Cable Type	Polyurethane outer
Connector Type	Thermocouple Plug
Connector Size	Miniature
Colour	Brown
Standards Met	Colour code (cores & jacket) to IEC-584-3
Thermocouple Type	T
Termination Type	Miniature Thermocouple Plug